



SMILE AS BRIGHT
AS THE SUNSHINE.

长沙奈斯医疗器械



CHANGSHA NICE MEDICAL INSTRUMENT CO.,LTD

COMPANY INTRODUCTION

Changsha Nice Medical Instrument Co., Ltd. is a professional orthodontic products manufacturer in Changsha of Hunan province, China.

It is a high-tech company specializing in the production, research and development, sales of orthodontic materials & products, dental treatment appliances, clinic tools and so on. Nice cooperated with China Central South University, and all developments are based on the strong dental clinical data of Hospital, the engineers and workers have more than 20 years experiences.

Nice also cooperated with the Laboratory of PMRI of Central South University (MIM technology pioneer in China, <https://pmri.csu.edu.cn/>) and the German N-Labor International Laboratory.

Nice is focus to offer the advance technology and design, high quality products, good service for every customers in the world.

Looking forward to build a long term business relationship with customers around the world, get to a win-win situation and a better future!

CONTENTS

Metal Self-ligating Bracket	01
Ceramic Self-ligating Bracket	07
Metal Bracket	09
Ceramic Bracket	13
Buccal Tube	15
Molar Band	17
Arch Wire	18
Orthodontic Accessories	20
Orthodontic Instruments	22
Clinic Burs	25

METAL SELF-LIGATING BRACKET—PERFECT



Low profile ,max improve patient comfortable.



Groove fillets are added under the four wings, which increases the space for the ligature wire or ligature ring ,and the same time, hold them well.



80 gauge vacuum diffusion encryption mesh makes the bonding stronger, laser marked in bottom.



Separate traction hook design, the doctor can choose whether insert the hook according to the needs.

With fashion triangle design for easy positioning and reliable fit.


Auxiliary Slot

Patented Locking structure, safe and reliable, easy to deal.

Large “R” round, reduce friction, improved the comfortable of patient.


Perfect Metal Self-ligating Bracket

Tooth Identification



How to open and close

To open, insert the open instrument into the slide release cavity and twist 90°.
To close, position the instrument on either side of the slot and squeeze the opposing prong.



Prescription	Hooks	Slot	Type
One Case Kit (Upper/Lower 5X5)	No Hooks	0.022	NC121-P1
	3 with Hooks	0.022	NC121-P2
	3, 4, 5 with Hooks	0.022	NC121-P3

Teech System		U1	U2	U3	U4	U5	L1	L2	L3	L4	L5
		L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R
MBT	Torque	+17	+10	0	-7	-7	-6	-6	-6	-12	-17
	Angulation	4	8	8	0	0	0	0	3	0	0
Roth	Torque	+12	+8	0	-7	-7	-2	-2	-11	-17	-22
	Angulation	5	9	11	0	0	0	0	5	0	0
Middle Torque	Torque	+15	+6	+7	-11	-11	-3	-3	+7	-12	-17
	Angulation	5	9	5	2	2	2	4	5	4	4
High Torque	Torque	+22	+13	+11	-11	-11	-3	-3	+13	-5	-17
	Angulation	5	9	5	2	2	2	4	5	4	4
Low Torque	Torque	+2	-5	-9	-11	-11	-11	-11	0	-12	-17
	Angulation	5	9	5	2	2	2	4	5	4	4

METAL SELF-LIGATING BRACKETS –POPULAR



Removable drop in smooth hooks



Auxiliary Slot



Laser welding reliable and secure



Color ID



Round edge and low profile design maximum patient comfort



Add one hook hole to adjust



Durable and strong clip easy open and close



Laser number mark on the bottom




METAL SELF-LIGATING BRACKET


Popular Metal Self-ligating Brackets

Tooth Identification

How to open and close



To open, insert the open instrument into the slide release cavity and twist 90°.
To close, position the instrument on either side of the slot and squeeze the opposing prong.



Prescription	Hooks	Slot	Type
One Case Kit	No Hooks	0.022/0.018	NC121-1
	3 with Hooks	0.022/0.018	NC121-2
	3, 4, 5 with Hooks	0.022/0.018	NC121-3

System \ Teeth		U1	U2	U3	U4	U5	L1	L2	L3	L4	L5
		L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R
MBT	Torque	+17	+10	-7	-7	-7	-6	-6	-6	-12	-17
	Angulation	4	8	8	0	0	0	0	3	0	0
Roth	Torque	+12	+8	0	-7	-7	-2	-2	-11	-17	-22
	Angulation	5	9	11	0	0	0	0	5	0	0
Middle Torque	Torque	+15	+6	+7	-11	-11	-3	-3	+7	-12	-17
	Angulation	5	9	5	2	2	2	4	5	4	4
High Torque	Torque	+22	+13	+11	-11	-11	-3	-3	+13	-5	-17
	Angulation	5	9	5	2	2	2	4	5	4	4
Low Torque	Torque	+2	-5	-9	-11	-11	-11	-11	0	-12	-17
	Angulation	5	9	5	2	2	2	4	5	4	4

METAL SELF-LIGATING BRACKETS-I

Removable drop in
smooth hooks

Laser welding
reliable and secure

Round edge and low
profile design maximum
patient comfort


Durable and strong clip
easy open and close

Laser number mark

METAL SELF-LIGATING BRACKETS-I


Metal Self-ligating Brackets

Tooth Identification



How to open and close

To open, insert the open instrument into the slide release cavity and twist 90°.
To close, position the instrument on either side of the slot and squeeze the opposing prong.



Prescription	Hooks	Slot	Type
One Case Kit	No Hooks	0.022	NC120-1
	3 with Hooks	0.022	NC120-2
	3, 4, 5 with Hooks	0.022	NC120-3

Teech System		U1	U2	U3	U4	U5	L1	L2	L3	L4	L5
		L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R
MBT	Torque	+17	+10	-7	-7	-7	-6	-6	-6	-12	-17
	Angulation	4	8	8	0	0	0	0	3	0	0
Roth	Torque	+12	+8	0	-7	-7	-2	-2	-11	-17	-22
	Angulation	5	9	11	0	0	0	0	5	0	0



CERAMIC SELF-LIGATING BRACKETS

- 1.High quality materials;
- 2.Thinner, smaller and more transparent;
- 3.Fine polishing,nice looking presentation;
- 4.Good design,preventing a loose pad or even a loose bracket from the teeth.

UR


UL

LR

LL


Nice Ceramic Self-ligating Bracket

Tooth Identification



How to open and close

To open, insert the open instrument into the slide release cavity and twist 90°.
To close, position the instrument on either side of the slot and squeeze the opposing prong.



Prescription	Hooks	Slot	Type
One Case Kit (Upper/Lower 5X5)	No Hooks	0.022	NC122-1
	3, 4, 5 with Hooks	0.022	NC122-2

System \ Teech		U1	U2	U3	U4	U5	L1	L2	L3	L4	L5
		L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R
MBT	Torque	+17	+10	-7	-7	-7	-6	-6	-6	-12	-17
	Angulation	4	8	8	0	0	0	0	3	2	2
Roth	Torque	+12	+8	0	-2	-7	-1	-1	-11	-17	-22
	Angulation	5	9	10	0	0	0	0	7	0	0




METAL BRACKETS

Produced with advanced Metal Injection Module (MIM) technique, Nice mini-metal bracket have smooth surface and low profile design, improve the comfort of patients and treatment efficiency.



Color Identification



Round wing tail curved design improves patient comfort and the arch wire is more stable.




Accurate slot data, easy to place the arch wire, and improve the efficiency of treatment.



Nice Mini Metal Brakcet

Tooth Identification



Prescription	Hooks	Slot	Type
One Case Kit (Upper/Lower 5X5)	No Hooks	0.022/0.018	NC0120-01
	3、4、5 with Hooks	0.022/0.018	NC0120-02


System \ Teech		U1	U2	U3	U4	U5	L1	L2	L3	L4	L5
		L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R
MBT	Torque	+17	+10	-7	-7	-7	-6	-6	-6	-12	-17
	Angulation	4	8	8	0	0	0	0	3	0	0
Roth	Torque	+12	+8	-2	-7	-7	-1	-1	-11	-17	-22
	Angulation	5	9	11	0	0	0	0	5	0	0
Edgewise	Torque	0	0	0	0	0	0	0	0	0	0
	Angulation	0	0	0	0	0	0	0	0	0	0

UR	UL
LR	LL

Smaller hook designed,ensure attachment of auxiliaries,providing maximum comfortable

Color Identification

Digital marking, easy to identify the tooth position.



Slot top hypotenuse structure design is convenient for placing the arch wire.


Slot top hypotenuse structure design is convenient for placing the arch wire.

METAL BRACKETS

Popular metal bracket use the top MIM technology,data is more accurate, the surface is smoother, the profile is smaller, the bonding is firmer and it's more comfortable for patient.


Popular Nice Metal Bracket

Tooth Identification



Prescription	Hooks	Slot	Type
One Case Kit (Upper/Lower 5X5)	No Hooks	0.022/0.018	NC0220-01
	3、4、5 with Hooks	0.022/0.018	NC0220-02

System \ Teech		U1	U2	U3	U4	U5	L1	L2	L3	L4	L5
		L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R
MBT	Torque	+17	+10	-7	-7	-7	-6	-6	-6	-12	-17
	Angulation	4	8	8	0	0	0	0	3	0	0
Roth	Torque	+12	+8	-2	-7	-7	-1	-1	-11	-17	-22
	Angulation	5	9	11	0	0	0	0	5	0	0
Edgewise	Torque	0	0	0	0	0	0	0	0	0	0
	Angulation	0	0	0	0	0	0	0	0	0	0




CERAMIC BRACKETS

- Use high MIM technology to promise the cute and smooth profile;
- The color is nearly the tooth, more pretty to wear;
- Accurate date to keep the treatment effectively;
- Two slot and three slot can be choosed.


Ceramic Bracket

Tooth Identification



Prescription	Hooks	Slot	Type
One Case Kit (Upper/Lower 5X5)	No Hooks	0.022/0.018	NC0320-01
	3, 4, 5 with Hooks	0.022/0.018	NC0320-02
Remarked: Two slot or three slot can be choosed.			

System \ Teech		U1	U2	U3	U4	U5	L1	L2	L3	L4	L5
		L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R	L/R
MBT	Torque	+17	+10	-7	-7	-7	-6	-6	-6	-12	-17
	Angulation	4	8	8	0	0	0	0	3	0	0
Roth	Torque	+12	+8	-2	-7	-7	-1	-1	-11	-17	-22
	Angulation	5	9	11	0	0	0	0	5	0	0
Edgewise	Torque	0	0	0	0	0	0	0	0	0	0
	Angulation	0	0	0	0	0	0	0	0	0	0



BUCCAL TUBE

- Wide entrance makes archwire insertion, easy to adjust archwire and reduces chair time;
- Top metal injection molding technology(MIM), ensures high accuracy data;
- Low profile design, smooth edge, improve the comfort of patients;
- Mesh compound contoured base provides a precise fit;
- Many styles and mesh shape can be choosed.





	System			Torq.	offsize	Item No.
Single Rectangle Buccal Tube Slot: 0.018/0.022"	Edgewise	Maxillary	1st Molar	0°	0°	NC-BT1061
			2nd Molar	0°	0°	NC-BT1071
		Mandibular	1st Molar	0°	0°	NC-BT1062
			2nd Molar	0°	0°	NC-BT1072
	ROTH	Maxillary	1st Molar	-14°	10°	NC-BT1063
			2nd Molar	-14°	10°	NC-BT1073
		Mandibular	1st Molar	-25°	4°	NC-BT1064
			2nd Molar	-25°	4°	NC-BT1074
	McBeTr Sliding	Maxillary	1st Molar	-14°	10°	NC-BT1065
			2nd Molar	-14°	10°	NC-BT1075
		Mandibular	1st Molar	-20°	0°	NC-BT1066
			2nd Molar	-20°	0°	NC-BT1076
Double Combination Buccal Tube (Slot and Round) Slot: 0.018/0.022" Round: Ø45"	Edgewise	Maxillary	1st Molar	0°	0°	NC-BT2061
		Mandibular	2nd Molar	0°	0°	NC-BT2062
	ROTH	Maxillary	1st Molar	-14°	10°	NC-BT2063
		Mandibular	2nd Molar	-30°	4°	NC-BT2064
	McBeTr Sliding	Maxillary	1st Molar	-14°	10°	NC-BT2065
		Mandibular	2nd Molar	20°	0°	NC-BT2066
Triple Combination Buccal Tube (Slot and Round) Slot: 0.018/0.022" Round: Ø45"	Edgewise	Maxillary	1st Molar	0°	0°	NC-BT3061
		Mandibular	2nd Molar	0°	0°	NC-BT3062
	ROTH	Maxillary	1st Molar	-14°	0°	NC-BT3063
		Mandibular	2nd Molar	-30°	0°	NC-BT3064
	McBeTr Sliding	Maxillary	1st Molar	-14°	0°	NC-BT3065
		Mandibular	2nd Molar	20°	0°	NC-BT3066

- High quality materials to promise the bands quality, laser markings for easy identification;
- Tempered with flattened interproximals for easier fitting;
- Broad buccal and lingual surfaces for greater retention;
- Different size and different combinations can be choosed;
- Pre-welded Buccal Tube use laser welding, make the entire structure strong and high efficiency.



First Molar Band


Spe.	30	30+	31	31+	32+	33	33+	34	34+	35	35+	36	36+	37	37+	38	38+	39	39+	40	40+	41	41+	42	42+	43
Perimeter (MM)	30.0	30.5	31.0	31.5	32.5	33.0	33.5	34.0	34.5	35.0	35.5	36.0	36.5	37.0	37.5	38.0	38.5	39.0	39.5	40.0	40.5	41.0	41.5	42.0	42.5	43.0


Second Molar Band

Spe.	S30	S30+	S31	S31+	S32+	S33	S33+	S34	S34+	S35	S35+	S36	S36+	S37	S37+	S38	S38+	S39	S39+	S40	S40+	S41	S41+	S42	S42+	S43
Perimeter (MM)	30.0	30.5	31.0	31.5	32.5	33.0	33.5	34.0	34.5	35.0	35.5	36.0	36.5	37.0	37.5	38.0	38.5	39.0	39.5	40.0	40.5	41.0	41.5	42.0	42.5	43.0

ARCH WIRE

- Nice Wire are manufactured by the high quality materials to promise the dimensional precision and surface finish;
- Continuous gentle force working for the whole treatment;
- Easy to treatment,it' s easy to bend into different shapes;
- Significant better results when working with nice bracket.





Item Name	Sep.	Size	Package
Niti Super Elastic Arch Wires	Round	0.012" ,0.014" ,0.016" ,0.018" ,0.020"	10 pcs/bag
	Rectangular	0.016" ×0.016" ,0.016" ×0.022" ; 0.017" ×0.022" ,0.017" ×0.025" ; 0.018" ×0.022" ,0.018" ×0.025" ; 0.019" ×0.025" ;0.021" ×0.025" .	10 pcs/bag
Niti Thermal Active Arch Wires	Round	0.012" ,0.014" ,0.016" ,0.018" ,0.020"	10 pcs/bag
	Rectangular	0.016" ×0.016" ,0.016" ×0.022" ; 0.017" ×0.022" ,0.017" ×0.025" ; 0.018" ×0.022" ,0.018" ×0.025" ; 0.019" ×0.025" ;0.021" ×0.025" .	10 pcs/bag
Stainless Steel Archwires	Round	0.012" ,0.014" ,0.016" ,0.018" ,0.020"	10 pcs/bag
	Rectangular	0.016" ×0.016" ,0.016" ×0.022" ; 0.017" ×0.022" ,0.017" ×0.025" ; 0.018" ×0.022" ,0.018" ×0.025" ; 0.019" ×0.025" ;0.021" ×0.025" .	10 pcs/bag
Copper NiTi Archwires	Round	0.012" ,0.014" ,0.016" ,0.018" ,0.020"	10 pcs/bag
	Rectangular	0.016" ×0.016" ,0.016" ×0.022" ; 0.017" ×0.022" ,0.017" ×0.025" ; 0.018" ×0.022" ,0.018" ×0.025" ; 0.019" ×0.025" ;0.021" ×0.025" .	10 pcs/bag
Niti Archwire Reverse curve	Round	0.012" ,0.014" ,0.016" ,0.018" ,0.020"	2 pcs/bag
	Rectangular	0.016" ×0.016" ,0.016" ×0.022" ; 0.017" ×0.022" ,0.017" ×0.025" ; 0.018" ×0.022" ,0.018" ×0.025" ; 0.019" ×0.025" ;0.021" ×0.025" .	2 pcs/bag

ORTHO ACCESSORIES (ELASTOMERICS)



Ligature Tie
1040 PCS/PACK 40 stricker*26 (round)

NC-E001



Cartoon Ligature Tie
1000 PCS/PACK Flower, Mickey and Kitty

NC-E002



Separator Elastomeric
15 PCS/PACK 26 O-ring in each pcs

NC-E003



Elastic Chain
15 FEETS/BOX, continue, short, middle long for choose

NC-E004



Orthodontic Elastics
5000 PCS/BOX,
Force: 2.5oz/3.5oz/4.5oz/6.5oz ;
D: 3.2mm/4.8mm/6.4mm/8mm/9.5mm

NC-E005



Ration Wedge
20 PCS/PACK


NC-E006



COLOR CHOOSD



ORTHO ACCESSORIES (ANOTHER)



Ligature Wire
NC-A001



Face Mask
NC-A002



Activity Crimpable Hook
NC-A003



Activity Archwire Stop
NC-A004



Crimpable Hooks
NC-A005



Crimpable Stop
NC-A006



Lingual Buttons
NC-A007



Button Chain
NC-A008



Open Coil Spring
NC-A009



Close Coil Spring
NC-A010



Dental Tray
NC-A011




Cheek Retractor
NC-A012



Orthodontic Mirror
NC-A013



Tooth Crown
NC-A014



ORTHODONTIC INSTRUMENTS





Light Wire Plier
NC-P101



Light Wire Cutter
NC-P102



Distal End Cutter
NC-P103



Torque Bending Plier
NC-P104



Arch Wire Cutter
NC-P105



Heavy Wire Cutter
NC-P106



KIM Plier
NC-P107



Tweed Plier
NC-P108



How's Plier
NC-P109



End Bending Plier
(Weingart Plier)
NC-P110



Bracket Removing Plier
NC-P111



Band Removing Plier
NC-P112



Adhesive Removing Plier
NC-P113



Crimpable Hooks Plier
NC-P114



Elastic Plier
NC-P115



Universal Plier
NC-P116

ORTHODONTIC INSTRUMENTS



Three cutter Plier
NC-P117



Vertical Notch Forming Plier
NC-P118



Ligature Tie Forming Plier
NC-P119



Omega Loop Bending Plier
NC-P120



Convertible Cap Removing Plier
NC-P121



Torque turn Plier
NC-P122



Mini Light Wire Plier
NC-P123



Long Weingart Plier
NC-P124



Three teeth Plier
NC-P125



Light Wire Forming Plier
NC-P126



Light Wire Forming Plier II
NC-P127



Lingual side Wire Forming Plier
NC-P128



Hook Fixing Plier
NC-P129



Double Tie Forming Plier
NC-130



Jarabak Plier
NC-131



End Bending Plier
NC-132

CLINIC BURS



- Use the most advanced MIM technology, and control the quality of the production process from materials to products strictly;
- All kinds of burs can be choosed;
- Top quality,easy to use and precise in treatment.

Clinic Burs

TAPER ROUND END CROSSCUT FISSURE SHAPE C33R

Rod Type	ISO No.	US No. Denburs No.	1700	1701	1702	Slot Length	Remarks
FG(19mm)	500 314 194 007	C33R.314	010	012	016	4.0/5.4	
FG xlg(25mm)	500 316 194 007	C33RL.316		012	016	8.0	

Repair Burs

ROUND SHAPE C1S

Rod Type	ISO No.	US No. Denburs No.	1/4S	1/2S	1S	2S	3S	4S	5S	6S	7S	8S	Slot Length	Cutter Number
FG(19mm)	500 314 001 006	C1S.314	005	006	008	010	012	014	016				\	12
RA(22mm)	500 204 001 006	C1S.204	005	006	008	010	012	014	016	018	021	023	\	12
RA lg(26mm)	500 205 001 006	C1S.205	005	006	008	010	012	014	016	018	021	023	\	12
FG xlg (25mm)	500 316 001 006	C1S.316	005	006	008	010	012	014	016				\	16

Repair Burs

STRAIGHT FLAT END FISSURE SHAPE C21S

Rod Type	ISO No.	US No. Denburs No.	57S	58S	59S	60S	61S	Slot Length	Cutter Number
FG xlg (25mm)	500 316 107 072	C21LN.316	010	012	014	016	018	8.0	12
FG xlg (25mm)	500 316 107 042	C21LM.316	010	012	014	016	018	8.0	20

Repair Burs

STRAIGHT ROUND END FISSURE SHAPE C21RS

Rod Type	ISO No.	US No. Denburs No.	1158S	1159S	1160S	1161S	Slot Length	Cutter Number
FG xlg (25mm)	500 316 137 006	C21RL.316	012	014	016	018	8.0	12
FG xxlg (28mm)	500 317 137 006	C21RL.317	012	014	016	018	8.0	12
FG xlg (25mm)	500 316 137 072	C21RLN.316	012	014	016	018	10.0	12
FG xlg (25mm)	500 316 137 042	C21RLM.316	012	014	016	018	10.0	20

Breaking Crown Burs

STRAIGHT ROUND END CROSSCUT FISSURE SHAPE C31LLR

Rod Type	ISO No.	US No. Denburs No.	1931	1932	Slot Length	Remarks
FG xlg (25mm)	500 316 141 008	C31LLRE.316	010	012	8.0	
FG xlg (25mm)	500 316 141 019	C31LLRX.316	010	012	8.0	

Breaking Crown Burs

TAPER ROUND END CROSSCUT FISSURE SHAPE C33LLR

Rod Type	ISO No.	US No. Denburs No.	1700	1701	1702	Slot Length	Remarks
FG xlg (25mm)	500 316 192 008	C33LLRE.316	012	014	016	10.0	
FG lg (21mm)	500 315 192 008	C33LRE.315	012	014	016	8.0	
FG xlg (25mm)	500 316 192 018	C33LRX.316	012	014	016	8.0	

Breaking Crown Burs

STRAIGHT ROUND END CROSSCUT FISSURE SHAPE C31RE

Rod Type	ISO No.	US No. Denburs No.	1931	1932	Slot Length	Remarks
FG(19mm)	500 314 137 008	C31RE.314	010	012	4.0	

Clinic Burs

TAPER FLAT END FISSURE SHAPE C23



Rod Type	ISO No.	US No. Denburs No.	168	169	169L	170	170L	171	171L	172	Slot Length	Remarks
FG (19mm)	500 314 168 006	C23.314	008	009	009	010	010	012	012	016	4.0/5.4	
FG xlg (25mm)	500 316 168 006	C23L.316				010		012			6.0	
HP (44.5mm)	500 104 168 006	C23L.104		009		010		012		016	6.0	

Clinic Burs

TAPER ROUND END FISSURE SHAPE C23R



Rod Type	ISO No.	US No. Denburs No.	1169	1169L	1170	1170L	1171	1171L	1172	Slot Length	Remarks
FG (19mm)	500 314 194 006	C23R.314	009	009	010	010	012	012	016	4.0/5.4	
RA (22mm)	500 204 194 006	C23R.204					012		016	4.0/5.4	

Clinic Burs

STRAIGHT FLAT END CROSSCUT FISSURE SHAPE C31



Rod Type	ISO No.	US No. Denburs No.	555	556	557	558	559	560	561	562	563	Slot Length	Remarks
FG (19mm)	500 314 107 007	C31.314		009	010	012	014	016				3.7/5.4	
FG xlg (25mm)	500 316 107 007	C31.316			010	012	014	016				5.4	
RA (22mm)	500 204 107 007	C31.204		009	010	012	014	016	018	031	023	3.7	
RA lg (26mm)	500 205 107 007	C31.205		009	010	012	014	016	018	031	023	5.4	
HP (44.5mm)	500 104 107 007	C31.104	008	009	010	012	014	016	018	021	023	6.0	

Clinic Burs

STRAIGHT FLAT END CROSSCUT FISSURE SHAPE C31L



Rod Type	ISO No.	US No. Denburs No.	556L	557L	558L	559L	560L	Slot Length	Remarks
FG (19mm)	500 314 110 007	C31L.314	008	010	012	014	016	5.5	
FG xlg (25mm)	500 316 110 007	C31L.316		010	012	014	016	8.0	
HP (44.5mm)	500 205 110 007	C31L.205		010	012			6.0	

Clinic Burs

STRAIGHT FLAT END CROSSCUT FISSURE SHAPE C31L



Rod Type	ISO No.	US No. Denburs No.	1555	1556	1557	1557L	1558	1558L	1559	1560	Slot Length	Remarks
FG (19mm)	500 314 137 007	C31R.314	008	009	010		012		014	016	3.7	
HP (44.5mm)	500 104 137 007	C31RL.104					012				6.0	
FG (19mm)	500 314 137 007	C31RL.314				010		012			5.5	
FG xlg (25mm)	500 316 137 007	C31RLL.316				010		012			10.0	

Clinic Burs

TAPER FLAT END CROSSCUT FISSURE SHAPE C33



Rod Type	ISO No.	US No. Denburs No.	699	699L	700	701	701L	702	702L	703	703L	Slot Length	Remarks
FG (19mm)	500 314 168 007	C33.314	009	009	010	012	012	016	016	021	021	4.0	
RA (22mm)	500 204 168 007	C33.204				012		016				5.4	
HP (44.5mm)	500 104 168 007	C33.104	009		010	012		016	016	021	021	5.4	

Clinic Burs

ROUND SHAPE C1



Rod Type	ISO No.	US No. Denburs No.	1/4	1/2	1	2	3	4	5	6	7	8	Slot Length	Remarks
FG (19mm)	500 314 001 001	C1.314	005	006	008	010	012	014	016				\	
FG xlg (25mm)	500 316 001 001	C1.316	005	006	008	010	012	014	016				\	
RA (22mm)	500 204 001 001	C1.204	005	006	008	010	012	014	016	018	021	023	\	
RA lg (26mm)	500 205 001 001	C1.205			008	010	012	014	016	018	021	023	\	
HP (44.5mm)	500 104 001 001	C1.104	005	006	008	010	012	014	016	018	021	023	\	

Clinic Burs

INVERTED CONE SHAPE C2



Rod Type	ISO No.	US No. Denburs No.	33 1/2	34	35	36	37	38	39			Slot Length	Remarks
FG (19mm)	500 314 010 001	C2.314	006	008	010	012	014	016				\	
RA (22mm)	500 204 010 001	C2.204	006	008	010	012	014	016	018			\	
HP (44.5mm)	500 104 010 001	C2.104		008	010	012	014	016	018	021	023	\	

Clinic Burs

STRAIGHT FLAT END FISSURE SHAPE C21



Rod Type	ISO No.	US No. Denburs No.	56	57	57L	58	59	38	61	62	63	Slot Length	Remarks
FG (19mm)	500 314 107 006	C21.314	008	010	010	012						4.0	
FG xlg (25mm)	500 316 107 006	C21.316		010		012	014	016				8.0	
RA (22mm)	500 204 107 006	C21.204		010		012	014	016	018			5.4	
HP (44.5mm)	500 104 107 006	C21.104	008	010		012	014	016	018	021	023	5.4	

Clinic Burs

PEAR SHAPE C7



Rod Type	ISO No.	US No. Denburs No.	329	330	331	58	332	332L	Slot Length	Remarks
FG (19mm)	500 314 237 001	C7.314	006	008	010		012		2.0	
FG (19mm)	500 314 238 001	C7L.314				010		012	4.0	
FG xlg (25mm)	500 316 237 001	C7.316		008	010				8.0	
RA (22mm)	500 204 237 001	C7.204			010				2.0	
RA lg (26mm)	500 205 237 001	C7.205						012	4.0	

Clinic Burs

STRAIGHT ROUND END FISSURE SHAPE C21R



Rod Type	ISO No.	US No. Denburs No.	1155	1156	1157	1158	1159	1160	1161	Slot Length	Remarks
FG (19mm)	500 314 137 006	C21R.314	008	009	010	012	014	016		4.0	
RA (22mm)	500 204 137 006	C21R.204	008		010	012		016	018	5.4	
RA lg (26mm)	500 205 137 006	C21R.205				012		016	018	6.0	

Repair Burs

TAPER ROUND END FISSURE SHAPE C23RS



Rod Type	ISO No.	US No. Denburs No.	1169S	1170S	1171S	1172S	Slot Length	Cutter Number
FG xlg (25mm)	500 316 194 072	C23RLLN.316	009	010	012	016	10.0	12
FG xlg (25mm)	500 316 194 042	C23RLLM.316	009	010	012	016	10.0	20
FG lg (21mm)	500 315 194 072	23RLLN.315	009	010	012	016	10.0	12

Surgical Tooth Extraction

ROUND SHAPE C1X



Rod Type	ISO No.	US No. Denburs No.	1X	2X	3X	4X	5X	6X	7X	8X	Slot Length	Remarks
FG (19mm)	500 314 001 019	C1X.314		010	012	014	016				\	
RA (22mm)	500 204 001 019	C1X.204	008	010	012	014	016	018	021	023	\	
RA lg (26mm)	500 205 001 019	C1X.205	008	010	012	014	016	018	021	023	\	

Surgical Tooth Extraction

TAPER FLAT END CROSSCUT FISSURE SHAPE C33L



Rod Type	ISO No.	US No. Denburs No.	700	701L	702L	Slot Length	Remarks
FG xlg (25mm)	500 316 171 007	C33L.316	010	012	016	6.0	
FG xxlg (28mm)	500 317 171 007	C33L.317	010	012	016	6.0	
FG xxxlg (32mm)	500 318 171 007	C33L.318	010	012	016	6.0	

TAPER FLAT END CROSSCUT FISSURE SHAPE C33F



Rod Type	ISO No.	US No. Denburs No.	1931	1932	Slot Length	Remarks
FG xlg (25mm)	500 316 172 298	C33LF.316	012	016	8.0	
FG xxlg (28mm)	500 317 170 298	C33F.317	012	016	4.0	
FG xxxlg (32mm)	500 318 171 298	33LF.318	012	016	6.0	

Minimally Invasive Burs

TAPER FLAT END FISSURE SHAPE C23LL



Rod Type	ISO No.	US No. Denburs No.	171	172	Slot Length	Cutter Number
FG xlg (25mm)	500 316 184 072	C23LLN.316	012	016	8.0	12
FG xlg (25mm)	500 316 184 042	C23LLM.316	012	016	8.0	20

Minimally Invasive Burs

TAPER FLAT END FISSURE SHAPE C23LL



Rod Type	ISO No.	US No. Denburs No.	171	172	Slot Length	Cutter Number
FG xlg (25mm)	500 316 184 071	C23LLZ.316	012	016	10.0	12
FG xlg (25mm)	500 316 184 074	C23LLK.316	012	016	10.0	12+1

	Ball drill (1)											
	US NO.		1/4	1/2	1	2	3	4	5	6	7	8
	Blade diameter (mm)		0.5	0.6	0.85	1.0	1.2	1.4	1.6	1.8	2.1	2.3
	ISO NO. 500 104 001 006		005	006	0085	010	012	014	016	018	021	023

	Inverted cone (2)									
	US NO.	33 ^{1/2}	34	35	36	37	38	39	40	41
	Blade diameter (mm)	0.6	0.8	1.0	1.2	1.4	1.6	1.8	2.1	2.3
	ISO NO. 500 104 010 006	006	008	010	012	014	016	018	021	023

	Cylindrical split groove (21)									
	US NO.		56	57	58	59	60	61	62	63
	Blade diameter (mm)		0.8	1.0	1.2	1.4	1.6	1.8	2.1	2.3
	ISO NO. 500 104 107 006		008	010	012	014	016	018	021	023

	Truncated cone split groove (23)	  		
	US NO.	169	170	171
	Blade diameter (mm)	0.9	1.0	1.2
	ISO NO. 500 104 168 006	009	010	012

	Cylindrical hemisphere split groove (21R)			
	US NO.	1156	1157	1158
	Blade diameter (mm)	0.9	1.0	1.2
	ISO NO. 500 104 137 006	009	010	012

	Cylindrical cross-filing (31)									
	US NO.	556	557	558	559	560	561	562	563	
	Blade diameter (mm)	0.9	1.0	1.2	1.4	1.6	1.8	2.1	2.3	
	ISO NO. 500 104 107 007	009	010	012	014	016	018	021	023	

	Crack drill (33)		A	A	A	A	A
	US NO.		699	700	701	702	703
	Blade diameter (mm)		0.9	1.0	1.2	1.6	2.1
	ISO NO. 500 104 168 007		009	010	012	016	021

	Wheel shape (3)	
	US NO.	14
	Blade diameter (mm)	1.2
	ISO NO. 500 104 040 006	0.12

PURSUE EXCELLENCE
ACHIEVE A SMILE





奈斯医疗
NICE MEDICAL

CHANGSHA NICE MEDICAL



长沙奈斯医疗
Changsha Nice Medical



Website: www.nicemedicalcs.com



E-mail: nicemedical@163.com



Tel: +86 15211074668 (whatsapp/facebook/linkedin)



Address: No. 501 in A9, Jinrongtongxin Industrial zone, Huizhi Road, Lugu Street, Yuelu District, Changsha City, Hunan Province, China.